

Medical Oxygen Cylinders

Medical Oxygen Cylinders must be inspected and pressure tested **every five years** by an accredited cylinder test house (Industrial Pressure Testing Ltd.)

Medical Oxygen has an **expiry period of three years**. Cylinder test periods and medical oxygen periods must be kept in line. Here are some examples.

Examples:

Cylinder retest date 2008/06 (June 2008)

The oxygen expiry date will be one year only as it must be kept in line with the next retest date of the cylinder.

Cylinder retest date 2010/06 (June 2010)

The oxygen expiry date will be three years (as normal)

Please check the Expiry Date of your Oxygen.(Batch Label) It may have expired.

Batch Label Example

BATCH: A14A9

DATE FILLED: 19/6/2007

EXPIRY DATE: 19/6/2010